



**TLSO  
BIVALVE**  
PDAC APPROVED  
• L0482



**LSO  
ANTERIOR  
OPENING  
& OVERLAP**  
PDAC APPROVED  
• L0640



**LSO  
BIVALVE  
SINGLE  
OPENING**  
PDAC APPROVED  
• L0640



**TLSO  
ANTERIOR  
OVERLAP  
& OPENING**  
PDAC APPROVED  
• L0486

**Available:** Custom-to-Measurement, Custom-to-Cast

(Orthosis modified by negative cast or Pre-modified Cast) Finished & Complete or Unfinished

**Anterior Opening:** Anterior Overlap or Anterior opening w/ tongue

**Bivalve Opening:** A over P, or P over A

**Overlap Options:** Smooth, Step (Butted) or with polymer tongues

**Anterior Options:** Flexion, neutral or pendulous

**Degree of Lordosis:** Custom

**Material Options:** LDPE, MPE, HDPE, Copolymer, Polypro, ABS, Others Available

**Liner Options:** Hypoallergenic Aliplast Foam 1/8", 3/16", 1/4", open cell breathable foam or unlined

**Closure Kit Attachments:** Single opening closure kit w/ 3 Straps, chafes, loops & Chicago Screws, Bivalve opening closure kit w/ 6 Straps, chafes, loops & Chicago Screws

**Sternal/Pectoral Pad Kit Attachments:** PPK (Pectoral Pad Kit) w/ straps, chafes, loops & screws, SPK (Sternal Pad Kit) w/ straps, chafes, loops & screws, Adjustable PPX (Pectoral Pad Kit) w/ straps, chafes, loops & screws, Adjustable SPX (Sternal Pad Kit) w/ straps, chafes, loops & screws

**Other Attachments:** Pneumatic lumbar bladder, Shoulder or axilla straps w/ chafes, loops & Chicago screws

**Custom Options:** Breast Buildup, Designer Pattern applied to brace & straps, Pubis Flare, Relief buildup, Trim Only, Waist pads, Other custom options available.

### Typical Indications

- Lumbar Fusion
- Degenerative, Herniated or Bulging Disc
- Lumbar Stabilization
- Post-Operative
- Spondylolisthesis
- Spinal Stenosis
- Wheelchair Support



MANUFACTURING SPECIALIST

Ph: 1-888-982-8181 • 1-770-513-7380 Fx: 1-877-288-4197 • 1-770-513-4365

[www.optecusa.com](http://www.optecusa.com)